



ASTIGMATISM


Binocular: to adjust the axis of the astigmatism of both eyes at the same time. The exam procedure is the same as the one with the classic targets (clock...) but without closing one of the two eyes.




Circular Optotypes: to adjust the power of the astigmatism using the Crossed Cylinders.




Clock: target for the identification of the astigmatism's axis.



Dial: target for the identification of the astigmatism's axis.



Dots: target for the identification of the astigmatism's axis.



Mono/Binocular: to adjust the astigmatism monocular power when both eyes are open. The exam procedure is the same as the one with the classic targets (clock. . .) but without closing one of the two eyes.

