



BALANCE


Balance: binocular polarized target to balance the near correction.




Balance for Kids: binocular polarized target with characters for kids to balance the near correction.



Binocular Radial: binocular target to correct the astigmatism. The exam procedure is the same as the one with the classic targets (clock...) but without closing one of the two eyes.



Osterberg Bichrom-Balance-Test: binocular test that unites the polarized filter to the red/green one, useful for the balance of the correction. Ask the patient what he/she sees, if the correction is balanced he/she should see the whole target without differences of resolution in the various parts of the target.



Red/Green: filter for the balance of the correction.

